



# Pergo®

Passion for Better



Pergo is a newly established medical company to design and manufacture innovative products based on more than 35 years of experience.

We allocate very large percentage of our turnover to our most important capital; human capital. We research, develop and manufacture patented single/multi-use, active and intelligent devices for global markets. We aim to be the most respected medical device design and manufacturing company in the World.

Pergo's vision is to be a followed as a world brand by adopting "creative", "fast", "reliable" and "irreproachable" principles to its all processes.



## Suction Liner System

- Highly flexible and transparent suction liner bag.
  - Super durability. Product would not damage even under toughest transport conditions or application.
  - Less transport cost. Approximately %50 less box dimension compared to the same quantity of traditional liners.
  - Less storage cost.
  - Less environmental waste.
  - Economic.
- Centralized and ergonomic removal handles.
  - Easy suction liner removal with one hand.
- When the suction liner is placed in the canister, it is ready to use. No other process is required to use the system.
- Fastest reach to 80 kPa vacuum in the market.
- Wide measuring scale for easy reading. Canister scale accuracy: +/- 100ml.
- Built-in brackets of canisters provides easy mounting on walls, rail supports and trolleys.
- Different suction tube options are available.
- FDA registered.
- CE marked.
- Latex free.

05TD0129	Suction Liner 2000 ml Canister Set
05TD0023	Suction Liner 2000 ml
05TD0061	Suction Liner 2000 ml with Solidifying Agent
05TD0013	Suction Liner Silicon Tube
05TD0126	Suction Liner Wall Hanger Set
05TD0128	Suction Liner Rail Hanger
05TD0136	Suction Liner Serial Connection Tube
05TD0177	Suction Liner Vacuum Control Valve EXT
05TD0176	Suction Liner Trolley - Stainless Steel
05TD0094	Suction Liner Extension Tube 180 cm with Connector
05TD0047	Suction Liner Vacuum T-Connector



## Skin Preparation Glove

- Easily collects shaved hair.
- Faster application with the large adhesive area.
- Does not irritate skin with the optimized adhesion.
- Can be applied both on flat and curved areas.
- Reduces surgical site contamination risk.
- Reduces post procedural infection risk.
- Avoids contamination with reverse removal design.
- Latex free.

05TD0003 Skin Preparation Glove



## Eye Safety Tape

- Easy, fast application and removal.
- Coloured, ergonomic and non-sticky handle.
- Designed to fit eye anatomy perfectly.
- Transparent structure allows 'easy pupil dilation control' at any time without removing the Eye Safety Tape.
- Does not hurt the skin, eyebrows and eyelashes while removing (post operation).
- Protects eyes from external splashes during the operation.
- Eyes do not dry out even during long operations.
- Medical grade material avoids any skin irritation.
- Reduces corneal abrasion risk.
- Hypoallergic, breathable and medical grade material.
- Patent protected.
- Latex free.

05TD0012 Eye Safety Tape Adult



## Clear Vision

- Light.
- Ergonomic.
- Antibacterial, reusable frame with different colors.
- Easily replaceable super clear disposable lens.
- Can be fitted on prescription eyewear and surgical masks.
- The large lens provides full protection for the eye and its surroundings.
- Anti-Fog.
- Latex free.
- Patent protected.

05TD0161	Clear Vision Lens DS (Single Use)
05TD0195	Clear Vision Frame Red (Reusable)
05TD0196	Clear Vision Frame Blue (Reusable)
05TD0198	Clear Vision Frame Black (Reusable)
05TD0200	Clear Vision Frame Yellow (Reusable)
05TD0202	Clear Vision Frame White (Reusable)
05TD0204	Clear Vision Frame Mixed Colours (Reusable)



## PostPartum Balloon

- 100% Silicone material allows for patient comfort.
- Pear Shaped balloon design is perfectly fit for uterine anatomy and shape.
- Super durable balloon.
- Does not require any surgical operations for removal.
- Dual lumen catheter allows inflating balloon and monitoring hemorrhage.
- Two syringes for inflating balloon rapidly.
- Bag Spike is equipped with check valve.
- Three way stopcock for easy, safe and rapid inflation.

05TD0206	PostPartum Balloon
----------	--------------------





## Mission

To develop, produce and distribute single and multiple use active and intelligent devices with an innovative understanding.

## Vision

Maintaining our presence as a world brand.

## Umbilical Catheter

- Catheter length: 30 cm.
- 3.5 Fr, 5 Fr and 7 Fr outer diameters.
- Can be used in Arterial and Ven catheterization.
- The raw material of the catheter is TPU (Thermosensitive Polyurethane).
- Thermosensitive material provides sufficient rigidity while entering the vein and prevents it from harming the veins by softening in the body.
- DEHP free.
- Markers spaced 1 cm apart between 4 cm and 25 cm.
- Rounded tip prevents irreparable damage to the vein.
- Radiopaque.
- Sterile and disposable.

05TD0225 Umbilical Catheter 3,5 Fr

05TD0268 Umbilical Catheter 5 Fr

05TD0269 Umbilical Catheter 7,5 Fr



## Probe Cover

- Reduces hospital infection rates.
- Offers additional sterile protection during procedures.
- Available in different lengths for complete coverage from transducer array to connector.
- Easy opening with special design for correct probe inserting.
- Guarding expensive imaging equipment during operation.
- Latex-Free material to protect latex-sensitive skins.

05TD0419 Sterile Probe Cover 14 x 25 cm

05TD0420 Sterile Probe Cover 14 x 40 cm

05TD0423 Sterile Probe Cover 14 x 80 cm

05TD0426 Sterile Probe Cover 14 x 120 cm

05TD0429 Sterile Probe Cover 14 x 160 cm

05TD0432 Sterile Probe Cover 14 x 250 cm



## Biopsy Valve

- Single use.
- Patented functional design.
- Compatible with Olympus, Pentax and Fujinon endoscopes.
- Reduces cross contamination.
- Leakage and splash minimization.
- Punched lid prevents back splash and ensures easy entry for accessories.
- Sterile / Non-sterile options available.
- FDA registered.
- Latex free.

05TD0010 Biopsy Valve Olympus-Fujinon Sterile

05TD0029 Biopsy Valve Pentax Sterile

05TD0054 Biopsy Valve Olympus-Fujinon

05TD0031 Biopsy Valve Pentax



## Endocage

- Patented Design.
- Semi-disposable.
- Re-openable lock if/when an accessory is left out by mistake.
- Easy hook onto the scope.
- Big inner volume.
- Product designed to allow maximum penetration of cleaning fluid.

05TD0130 Endocage





## ERCP Drape

- Four pockets for forceps and tools.
- Fixed from 3 different location for high stability.
- High liquid absorption with a special structure.
- Extra strong adhesive fixation tapes.
- Latex free.
- Sterile.

05TD0205 ERCP Drape



## Co2 Tubing Set

- Universal design.
- Works with any disposable water bottle system.
- Works effectively with any insufflator on the market.
- Works with Olympus, Fujinon and Pentax endoscopes.
- Simple Luer connections.
- Tubing can be disposed of at the end of each day.
- Features a hydrophobic filter to prevent back-flow of water into the insufflator.

05TD0165 Endoscopy Insufflation CO2 Tubing Set





Pergo Medikal ve İlaç Sanayi A.Ş.  
Buca OSB Mah. 2/20 Sok. BEGOS Sit. No:9  
Buca - İzmir / Turkey  
☎ +90 232 2908929 ✉ cservices@pergomed.com  
🌐 www.pergomed.com