

CLEAR VISION COMFORT FACE SHIELD




Intended Use

Clear Vision Comfort Face Shield is designed and developed to protect face from blood, other body fluid splashes during medical interventions or practices. It also provides a clear vision to the user. Product is used to prevent contaminations through eyes, nose, mouth and to reduce the risk of diseases from spreading.

Key Features

- Light,
- Ergonomic,
- Antifog,
- Single use,
- Can be used along with eyeglasses and masks,
- Wide lens provides full face protection,
- Elastic band provides a stable position on head,
- Soft forehead sponge provides long hours use of the face mask without any discomfort,
- Double sided protective layers prevent the film against scratches and dirt, can be easily peeled of before use.
- Face Shield is presented in four versions with Standart PET/PC Lens and Long PET / PC Lens options,

System Parts



Models

REF. No	Product Name	Barcode No	Description
05TD3054	Clear Vision Comfort Face Shield	8681756018666	PC lens, double side protective film, single use.
05TD3070	Clear Vision Comfort Face Shield L	8681756018727	PC long lens, double side protective film, single use.
05TD3119	Clear Vision Comfort Face Shiled Eco	8681756018857	PET lens, double side protective film, single use.
05TD3120	Clear Vision Comfort Face Shiled Eco L	8681756018871	PET long lens, double side protective film, single use.



Precautions

- Do not detach the parts of the Face Shield.
- Store the Face Shield in a dry environment.
- Face Shield is single use. Do not reuse.
- Do not clean the Face Shield by any detergent or chemicals. Dispose properly after use.

Quality

- Manufactured under ISO 13485:2016 QMS.

Shelf Life

- 5 Years.

