

# Co2 TUBING SET

---

ENDOSCOPY



## Intended Use

Adequate distension of the GI lumen is required for safe advancement of endoscopes and for careful visualisation of the mucosa. CO<sub>2</sub> is a superior alternative to room air for insufflation during GI Endoscopy. CO<sub>2</sub> pressure and flow rate is controlled by insufflators and the gas is passed through a water bottle that is connected to the endoscope by tubing.

CO<sub>2</sub> Tubing Set is used to connect a CO<sub>2</sub> insufflator to the water bottle of the endoscope.



## Key Features

- Universal design.
- Works with any disposable water bottle system.
- Works effectively with any insufflator on the market.
- Works with Olympus, Fujinon and Pentax endoscopes.
- Simple Luer connections.
- Tubing can be disposed of at the end of each day.
- Features a hydrophobic filter to prevent back-flow of water into the insufflator.

## Models

REF 05TD0165 Endoscopy Insufflation CO<sub>2</sub> Tubing Set

Pergo Medikal ve İlaç Sanayi A.Ş.  
Buca OSB Mah. 2/20 Sok. BEGOS Sit. No:9  
Buca - İzmir / Turkey

☎ +90 232 2908929 ✉ cservices@pergomed.com  
🌐 www.pergomed.com