

ENDOCAGE

ENDOSCOPY



Intended Use

Endocage, is a semi-disposable container used for safe transport and traceability of endoscopy accessories during disinfection and storage. Product guarantees that accessories are aligned with the relevant endoscope. Product is designed to fit four accessories.

Key Features

- Patented Design.
- Semi-disposable.
- Re-openable lock if/when an accessory is left out by mistake.
- Easy hook onto the scope.
- Big inner volume.
- Product designed to allow maximum penetration of cleaning fluid.

Models

REF 05TD0130 Endocage



Pergo Medikal ve İlaç Sanayi A.Ş.
Buca OSB Mah. 2/20 Sok. BEGOS Sit. No:9
Buca - İzmir / Turkey

☎ +90 232 2908929 ✉ cservices@pergomed.com
🌐 www.pergomed.com