


# EZEUSE (ENDOSCOPE TRAY BAG)

---

## ENDOSCOPY



## Intended Use

- Transporting the sterile endoscope to the operating room after endoscope sterilization.
- Transporting the contaminated endoscope to the endoscope sterilization after endoscopy procedure at the operating room.
- Sterile endoscope storage and transport bag.

## Components of The Product

- Green and Red square bag with two adhesive strips,
- Inner Liner.

## Product Features for Ezeuse Bag

- Universal bag for every carrying tray.
- Special sterile endoscope storage and reversible carrying bag (and liner).
- Transparent window to indentify contents easily.
- Clean/contaminated markings to easily distinguish as clean or dirty.
- Completely contains the soiled equipment, reducing the risk of spreading infection.
- Sealing of the bag prevents contamination.



REF NO: 05TD0812

## Quality

- Manufactured under ISO 13485 Quality Management System
- CE marked according to 93/42/EEC directive and classified as Class I sterile medical device.

## Biocompatibility

- The product does not incorporate substances of human blood derivative, animal originated tissues, phthalates, medicinal product, latex, radioactive material and electromagnetic waves.
- The product has no contact with patient.

## Sterilization

- Ezeuse is sterilized with Ethylene oxide.

## Storage Conditions

- Must be stored in dry, ventilated environments, away from direct sunlight.

## Shelf Life

- 5 years

Pergo Medikal ve İlaç Sanayi A.Ş.  
Buca OSB Mah. 2/20 Sok. BEGOS Sit. No:9  
Buca - İzmir / Turkey

☎ +90 232 2908929 ✉ cservices@pergomed.com  
🌐 www.pergomed.com