

POSTPARTUM BALLOON

GYNECOLOGY



Intended Use

Postpartum Balloon is a disposable, multiple lumen catheter attached to an inflatable balloon system designed to provide tamponade for controlling hemorrhage from the uterus and vagina.

Features

- 100% Silicone material allows for patient comfort.
- Pear Shaped balloon design is perfectly fit for uterine anatomy and shape.
- Super durable balloon.
- Does not require any surgical operations for removal.
- Dual lumen catheter allows inflating balloon and monitoring hemorrhage.
- Two syringes for inflating balloon rapidly.
- Bag Spike is equipped with check valve.
- Three way stopcock for easy, safe and rapid inflation.
- Balloon can intake upto 5000 ml saline water. Considerably superior than other available products in the market.



Models

REF 05TD0206 Postpartum Balloon

Pergo Medikal ve İlaç Sanayi A.Ş.
Buca OSB Mah. 2/20 Sok. BEGOS Sit. No:9
Buca - İzmir / Turkey

☎ +90 232 2908929 ✉ cservices@pergomed.com
🌐 www.pergomed.com