



# Passion for Better



**Pergo** is a dynamic medical technology company established on **more than 35 years** industry experience.

We allocate large percentage of our income to our most important capital; human capital.

We research, develop and manufacture patented single/multi-use, active and intelligent devices for global markets. We aim to be the most respected medical device design and manufacturing company in the World.

## **Pergo Principles**

- Creative
  - **Fast**
- Irreproachable
  - **Reliable**

# INDEX

4 Suction Management


6 Operating Room

9 Surgical Equipment Covers

15 Endoscopy

19 Gynecology

19 Breathing Therapy





05TD0013



05TD0136




05TD0126



05TD0128



05TD0454



05TD1974



05TD0047



05TD0176



05TD0405



05TD0854

## Suction Liner System

- 100% inline robotic quality controlled production. No defected product leaves the factory.
- Highly flexible and transparent suction liner bag,
  - Super durability. Product would not damage even under toughest transport conditions or application.
  - Less transport cost. Approximately 65% less box dimension compared to the same quantity of traditional liners.
  - Less storage cost.
  - Less waste to the environment.
  - Economic.
- Centralized and ergonomic removal handles,
  - Easy suction liner removal with only one hand.
- As soon as suction liner is assembled onto the canister, system is ready to use. No additional vacuum connection required.
- Wide measuring scale for easy reading.
- Built-in brackets of canisters provides easy mounting on walls, rail supports and trolleys.
- Patented folding system,
  - Easy set up of liner into canister.
  - Easy removal of liner from the canister.
- FDA listed.
- CE marked.
- Latex free.

**05TD0852** Suction Liner HF 1000 ml (Sample Port)  
**05TD1665** Suction Liner HF 1500 ml (Sample Port)  
**05TD0405** Suction Liner HF 2000 ml (Sample Port)  
**05TD0939** Suction Liner HF 3000 ml (Sample Port)  
**05TD0853** Suction Liner Canister HF 1000 ml  
**05TD1506** Suction Liner Canister HF 1500 ml  
**05TD0854** Suction Liner Canister HF 2000 ml  
**05TD0855** Suction Liner Canister HF 3000 ml

**05TD0013** Suction Liner Silicone Tube  
**05TD0126** Suction Liner Wall Hanger Set  
**05TD0128** Suction Liner Rail Hanger  
**05TD0136** Suction Liner Serial Connection Tube  
**05TD0454** Suction Liner Vacuum Control Valve L  
**05TD0176** Suction Liner Trolley - Stainless Steel  
**05TD1974** Suction Liner Extension Tube 180 cm (Female connector, Kapkon connector)  
**05TD0047** Suction Liner Vacuum T Connector



Please refer to the brochures for all other codes related to the products.

## Eye Safety Type


- Easy, fast application and removal.
- Coloured, ergonomic and non-sticky handle.
- Designed to fit eye anatomy perfectly.
- Transparent structure allows 'easy pupil dilation control' at any time without removing the Eye Safety Tape.
- Does not hurt the skin, eyebrows and eyelashes while removing (post operation).
- Protects eyes from external splashes during the operation.
- Eyes do not dry out even during long (4 hours +) operations.
- Medical grade material avoids any skin irritation.
- Reduces corneal abrasion risk.
- Hypoallergenic, breathable and medical grade material.
- Patent protected.
- Latex free.



- 05TD0012** Eye Safety Tape Adult
- 05TD0020** Eye Safety Tape Pediatric
- 05TD1196** Eye Safety Tape Dispenser

## Eye Safety Tape Adult Sensitive

- Designed for patients with sensitive skin.
- Suitable for skin prone to developing adhesive allergies.
- Suitable for elderly patients with delicate skin.
- Reduced risk of allergy.
- Waterproof tape structure.
- Latex free.



- 05TD0019** Eye Safety Tape Adult Sensitive

## Skin Preparation Glove

- Easily collects shaved hair.
- Faster application with the large adhesive area.
- Does not irritate skin with the optimized medical adhesion.
- Can be applied both on flat and curved areas.
- Reduces surgical site contamination risk.
- Reduces post procedural infection risk.
- Avoids contamination with reverse removal design.
- Latex free.

**05TD0003** Skin Preparation Glove




## Clear Vision

- Light.
- Ergonomic.
- Antibacterial, reusable frame with different colors.
- Easily replaceable super clear disposable lens.
- Can be fitted onto prescription eyewear and surgical masks.
- Large lens provides full protection for the eye and its surroundings.
- Anti-Fog.
- Patent protected.
- Latex free.

**05TD0161** Clear Vision Lens DS (Single Use)

**05TD0204** Clear Vision Frame Mixed Colours (Reusable)





## Probe Cover

- Reduces hospital infection rates.
- Offers additional sterile protection during procedures.
- Available in different lengths for complete coverage from transducer array to connector.
- Easy application thanks to telescopic folding.
- Patented Adhesive models are used without ultrasonic gel between the transducer and the drape. No interference with ultrasonic image.
- Telescopic folded by robots without any human contact. Lowest possible bioburden level before sterilization.
- Probe covers in different lengths and kits including 20ml. sterile gel are available upon request.
- Sterile.
- Latex free.

**05TD0083** Probe Cover with Adhesive 18x80 cm  
**05TD0676** Probe Cover with Adhesive 18x120 cm  
**05TD0678** Probe Cover with Adhesive 18x160 cm

**05TD0423** Probe Cover 14x80 cm  
**05TD0426** Probe Cover 14x120 cm  
**05TD0429** Probe Cover 14x160 cm



Please refer to the brochures for all other codes related to the products.

## Tee Probe Cover

- Compatible with all standard TEE probes.
- Provides full coverage from the probe tip to be connector.
- Sterile.
- Latex free.

**05TD7190** Tee Probe Cover 3,5x122 cm Set  
**05TD7196** Tee Probe Cover 3,5x122 cm



## Mini C-Arm Cover

- It completely surrounds the equipment and prevents contamination.
- Special U-shape design conforms easily to the device.
- Three durable straps to properly secure the device to the machine.
- Compatible with a variety of mini C-arm models.
- Versions without a rubber band are available upon request.
- Sterile.
- Latex free.

**05TD7003** Mini C-Arm with Rubber Band 138x160 cm

**05TD8972** Mini C-Arm with Rubber Band 138x220 cm



## Standard Camera Cover

- Provides wide device compatibility with different length and width options.
- Safely isolates the camera system.
- Easy to use and apply with a transparent and triangle perforated tip design.
- Easy application with central protuberance and adhesive strip.
- Sterile.
- Latex free.

**05TD0365** Telescopic Folded Camera Cover 13x250 cm



## Cardboard Camera Cover

- Specially designed triangle perforated tip provides easy application.
- Its telescopically folded structure makes it easier and faster to cover the camera.
- Compatible with all camera models.
- It can be produced in different sizes upon request.
- Provides wide device compatibility with different length and width options.
- Safely isolates the camera system.
- Easy application with cardboard applicator with hole and adhesive strip.
- Sterile.
- Latex free.

**05TD1294** Cardboard Camera Cover 13x250 cm

**05TD0135** Cardboard Camera Cover 18x250 cm



## Universal Light Handle Cover

- Mostly compatible with different operating room ceiling lights of various sizes (from Ø20 mm to Ø40 mm).
- Secure locking tabs keep the cover firmly in place, preventing any risk of slipping or falling during use.
- Easy push fit mounting with one hand.
- Standard (140 mm) and Long versions (165 mm) are available.
- Sterile.
- Latex free.

**05TD0032** Universal Light Handle Cover

**05TD9544** Universal Light Handle Cover (Long)



## Ring Camera Cover

- Specially folded for fast application.
- Provides wide device compatibility.
- Safely isolates the camera system.
- Easy to use and apply with a transparent and triangle perforated tip design.
- Easy application with ring applicator and adhesive strip.
- Sterile.
- Latex free.

**05TD3573** Ring Camera Cover 14x250 cm



## Camera Cover with Elastic Tip

- Specially designed elastic tip provides easy application.
- Its telescopically folded structure makes it easier and faster to cover the camera.
- Compatible with all endoscopic camera models.
- It can be produced in different sizes upon request.
- Sterile.
- Latex free.

**05TD9243** Camera Cover 15x300 cm (Elastic Tip)



## C-Arm Cover (With Clips)

- Reduces hospital infection rates.
- Includes two scopy drapes and one C Arm Drape.
- Models including scopy covers in different sizes and quantities are available upon request.
- Scopy arm is easily covered by the clip-on cover and prevents contamination.
- Easy to put up and take down thanks to the flexible clips.
- Perfect surroundings to equipment NO contamination.
- Compatible with Siemens Philips, Ziehm and GE C-Arm machines.
- Sterile.
- Latex free.

**05TD1218** C-Arm Cover Large

**05TD1989** C-Arm Cover Medium

**05TD0618** C-Arm Cover Small



## C-Arm Cover (Without Clips)

- Specially folded for easy application.
- It consists of 1 unit C arm cover.
- It can be produced in different sizes upon request.
- Sterile.
- Latex free.

**05TD3318** C-Arm Long Cover 90x225 cm



## Dome Bag

- Compatibility with different devices can be achieved thanks to the variety of sizes.
- It can be produced in different sizes upon request.
- Sterile.
- Latex free.

- 05TD9068** Dome Bag 130x65 cm
- 05TD9478** Dome Bag 85x42,5 cm
- 05TD9479** Dome Bag 95x47,5 cm
- 05TD9481** Dome Bag 140x70 cm



## Banded Bag

- Compatibility with different devices can be achieved thanks to the variety of sizes.
- It can be produced in different sizes upon request.
- Double sided adhesive tapes attached to collect and secure excess cover material.
- Sterile.
- Latex free.

- 05TD3679** Banded Bag 75x75 cm
- 05TD3680** Banded Bag 75x90 cm
- 05TD3681** Banded Bag 90x50 cm
- 05TD1432** Banded Bag 90x75 cm
- 05TD9418** Banded Bag 80x100 cm
- 05TD9419** Banded Bag 60x75 cm






## Biopsy Valve

- Single use.
- Patented functional design.
- Compatible with Olympus and Fujinon endoscopes.
- Reduces cross contamination risk.
- Leakage and splash minimization.
- Sterile / Non-sterile options available.
- Latex free.

- 05TD0010** Biopsy Valve Olympus-Fujinon Sterile
- 05TD0054** Biopsy Valve Olympus-Fujinon
- 05TD0364** Biopsy Valve Olympus-Fujinon Non-Steril Single Package



## Bite Block

- Single use.
- Strap to secure comfortably in place during the procedure.
- Compatible with mouth anatomy.
- Fits different sizes of endoscopes.
- Non-sterile.
- Latex free.

- 05TD0777** Endoscopy Adult Bite Block with Strap



## Endoscope Valve Set

- Patent pending for Air-Water and Suction Valves , Patented design for Biopsy Valve
- Compatible with 140/160/180/190/240/260 Olympus GI Endoscopes.
- Highest suction rate among all other single use Valve Sets.
- Color coded Air-Water and Suction Valves.
- Easy removal of Air/Water and Suction Valves from the scope via outer ring form.
- Superior durability of Air/Water and Suction Valves more than 1000 times pressing.
- Easy mounting and dismounting of Water Jet Connector to the scope via optimized design.
- Perfect backflow prevention performance up to 4 bar negative pressure for Water Jet Connector.
- Latex free.

- 05TD1191** Air/Water and Suction Valve Sterile (2-piece kit)
- 05TD3607** Air/Water, Suction and Biopsy Valve Sterile (3-piece kit)
- 05TD1636** Air/Water, Suction, Biopsy Valve and Auxiliary Water Jet Connector Sterile (4-piece kit)



## Polyp Trap

- Designed to facilitate uninterrupted suction and provide separation of polyps from collected endoscopic effluent.
- Clear plastic construction allows for quick visual verification of polyp specimens.
- Single use to prevent cross-contamination.
- Non-sterile.
- Latex free.



**05TD0803** P-Trap

## Ezeuse Endoscope Tray Bag

- Universal bag for every carrying tray.
- Transparent window for easy identification of contents.
- Clean/contaminated markings to easily distinguish as clean or dirty.
- Completely contains soiled equipment, reducing the risk of spreading infection.
- Sealing of the bag prevents contamination.
- Latex free.



**05TD0812** Ezeuse

## Co2 Tubing Set


- Universal design.
- Compatible with any disposable water bottle system.
- Works effectively with any insufflator in the market.
- Compatible with Olympus, Fujinon and Pentax endoscopes.
- Simple set up luer connections.
- Tubing can be disposed at the end of each day.
- Features a hydrophobic filter to prevent back-flow of water into the insufflator.
- Latex free.



**05TD0165** Endoscopy Insufflation CO2 Tubing Set

## PostPartum Balloon

- 100% Silicone balloon and catheter for maximum patient comfort.
- Pear shaped balloon perfectly fits to uterine anatomy.
- Dual lumen catheter allows inflating balloon and monitoring hemorrhage.
- Set includes two syringes for rapid inflation of balloon.
- Bag spike is equipped with check valve.
- Custom design three way stopcock for easy, safe and rapid inflation of balloon.
- Maximum inflation volume is 500 ml.
- Sterile.
- Latex free.



05TD0206 PostPartum Balloon

## BREATHING THERAPY

### Aergo Breathing Therapy Devices

- Supports breathing with vPEP technology.
- Helps clear secretions and keep airways open.
- Adjustable resistance fits different patient needs.
- Portable, ergonomic, and easy to clean.
- 100% QC tested for stable airflow and consistent performance.
- Enables nebulizer use with the Aergo Pro option.



05TD10963 Aergo One Breathing Therapy Device



05TD10964 Aergo Pro Breathing Therapy Device



Pergo Medikal ve İlaç Sanayi A.Ş.  
Buca OSB Mah. 2/20 Sok. BEGOS Sit. No:9  
Buca - İzmir / Turkey

+90 232 2908929  [cservices@pergomed.com](mailto:cservices@pergomed.com)

 [www.pergomed.com](http://www.pergomed.com)