

SUCTION LINER SYSTEM

SUCTION MANAGEMENT



SUCTION MANAGEMENT

Key Features:

- Highly flexible and transparent suction liner bag.
 - Super durability. Product would not damage even under toughest transport conditions or application.
 - Less transport cost. Approximately %50 less box dimension compared to the same quantity of traditional liners.
 - Less storage cost.
 - Less environmental waste.
 - Economic.
- Centralized and ergonomic removal handles.
 - Easy suction liner removal with one hand.
- No additional vacuum connection is required. As soon as liner placed into the canister, system is ready to use.
- Lowest system set up time among all suction liners. 5 seconds (Pergo) vs 18.5 seconds (Hospira, Cardinal)
- Fastest reach to 80 kPa vacuum in the market.
- Wide measuring scale for easy reading.
Canister scale accuracy: +/- 100ml.
- Built-in brackets of canisters provides easy mounting on walls, rail supports and trolleys.
- Different suction tube options are available.
- Canister design according to ISO 19054:2005
- FDA registered.
- CE marked.
- Latex free.



SUCTION MANAGEMENT

System Parts

REF	05TD0129	Suction Liner 2000 ml Canister Set
REF	05TD0023	Suction Liner 2000 ml
REF	05TD0061	Suction Liner 2000 ml with Solidifying Agent
REF	05TD0013	Suction Liner Silicon Tube
REF	05TD0126	Suction Liner Wall Hanger Set
REF	05TD0128	Suction Liner Rail Hanger
REF	05TD0136	Suction Liner Serial Connection Tube
REF	05TD0454	Suction Liner Vacuum Control Valve L
REF	05TD0176	Suction Liner Trolley - Stainless Steel
REF	05TD0094	Suction Liner Extension Tube 180 cm with Connector
REF	05TD0047	Suction Liner Vacuum T-Connector



05TD0129



05TD0023



05TD0094



05TD0013



05TD0126



05TD0128



05TD0136



05TD0047



05TD0454



05TD0176

Pergo Medikal ve İlaç Sanayi A.Ş.
Buca OSB Mah. 2/20 Sok. BEGOS Sit. No:9
Buca - İzmir / Turkey